



FSU Biohazardous

Non-Infectious Waste

Non-pathogenic bacteria



Place in **clear, non-labeled autoclave bag**



Waste that is not ready to be autoclaved, ensure that it is in a closed container labeled **"General waste to be autoclaved only"**



Place the **clear bag** in an **autoclaveable container** and autoclave on appropriate settings

Infectious Waste

Disease causing organisms, human pathogen containing waste



Place a **red labeled biohazard autoclave bag**



Autoclave on appropriate settings

Sharps

Broken/Unbroken Glass, Pipette Tips, Pasteur Pipettes

Razor Blades, Syringes, Scalpels, Needles



Uncontaminated place in cardboard box labeled **"Broken Glass"**



Contaminated, place in puncture resistant labeled container



Place in **red puncture resistant container**



Can be discarded as **Trash**



Put in **Biomedical Waste Box**, discard as **Biomedical Waste** (Not to exceed 30 days)

Contact Richard Le or Janice Dodge